

**RAKIRO BIOTECH SYSTEMS PVT LTD**

R-466, TTC Industrial Area, MIDC  
Rabale, Navi Mumbai - 400701  
Tel No. +91-022-47804040  
Email :- sales@rakiro.net

**Doc No :** TDSAIPPH04  
**Date :** 01-02-2024  
**Type :** AQUASOL INDICATOR PAPER  
**Product Code:** AIPPH04

**PRODUCT DATA SHEET****1 INFORMATION**

Product Code: AIPPH04  
Product Name: AQUASOL - INDICATOR PAPER NARROW RANGE pH 3.5 6.0 (10 X 20 STRIPS)

**2 DESCRIPTION**

Ph Indicator Paper Testing have been developed with a mixture of several indicators, which on contact with the sample solution assume a particular color. A check against the color comparison table supplied allows the parameter to be determined. AQUASOL range offers a range of pH indicator and other test papers to meet your specific testing needs. Our AQUASOL Indicator Paper combine ease-of-use with accuracy, consistency and unmatched quality. The convenience of using AQUASOL Indicator Paper for the rapid determination of pH and other parameters has many applications in laboratories and industry.

**3 SALIENT FEATURES**

- Instant readings
- Economical in Kangaroo Pack
- Accurate for a wide range of testing
- Convenient and portable

**4 ANALYTICAL DATA**

Appearance : 10 Booklet of 20 Paper Strip each in Plastic Container  
Parameter : pH  
Range : 3.5 - 6.0

**5 HANDLING INSTRUCTIONS**

Keep away from Direct Heat, Sunlight and Moisture. Keep the container closed when not in use. The product should not be swallowed and prolonged contact with the skin should be avoided.

**6 Direction to use :**

- 1) Tear off a strip from the booklet and insert into the solution to be tested.
- 2) For highly viscous or stained liquids, the substance is dripped onto the paper.
- 3) Compare with the colour chart.

**Note :**

- 1) No Liability accepted for the accidents arising in handling or use.
- 2) Keep away from heat and sunlight.
- 3) Store in cool dry place & protect from chemical fumes.
- 4) Any discolouration in dry new paper does not impair the efficacy of the product.